

# GABRG1 antibody

Catalog Number: 110812

## Product name

GABRG1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GABRG1 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive WB detected in mouse eye tissue. Observed molecular weight by Western-blot: 54 kDa

## Preparation

This antibody was obtained by immunization of GABRG1 recombinant protein (Accession Number: NM\_173536). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:110812(GABRG1 antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:110812(GABRG1 antibody) at dilution of 1:100 (under 40x lens)