

# GAD1 antibody



Catalog Number: 110816

## Product name

GAD1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GAD1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, mouse brain tissue. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 67 kDa, 25 kDa

## Preparation

This antibody was obtained by immunization of GAD1 recombinant protein (Accession Number: XM\_017003758). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110816(GAD1 Antibody) at dilution of 1:600



IP Result of anti-GAD1 (IP:Catalog No:110816, 4ug; Detection:Catalog No:110816 1:500) with HEK-293 cells lysate 3680ug.