

# GADD45G antibody



Catalog Number: 110821

## Product name

GADD45G antibody

## Specificity

Human; other species not tested.

## Antibody description

GADD45G Rabbit Polyclonal antibody. Positive WB detected in UV treated HeLa cells, human placenta tissue. Observed molecular weight by Western-blot: 17kd

## Preparation

This antibody was obtained by immunization of GADD45G recombinant protein (Accession Number: NM\_006705). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



UV treated HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110821(GADD45G Antibody) at dilution of 1:300