

Galc antibody

Catalog Number: 110832

Product name

Galc antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Galc Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, A375 cells. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 80kd, 30kd

Preparation

This antibody was obtained by immunization of Galc recombinant protein (Accession Number: BC086671). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110832(Galc antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:110832(Galc antibody) at dilution of 1:50 (under 10x lens)