

# Galectin 2 antibody



Catalog Number: 110836

## Product name

Galectin 2 antibody

## Immunogen

[Human LGALS2 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

Galectin 2 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in human colon tissue. Observed molecular weight by Western-blot: 15kd

## Preparation

This antibody was obtained by immunization of Galectin 2 recombinant protein (Accession Number: BC029063). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



human colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:110836(LGALS2 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:110836(LGALS2 antibody) at dilution of 1:100 (under 10x lens)