

Galectin-7 antibody



Catalog Number: 110841

Product name

Galectin-7 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Galectin-7 Rabbit Polyclonal antibody. Positive WB detected in mouse skin tissue, mouse small intestine tissue. Positive IHC detected in human cervix tissue. Observed molecular weight by Western-blot: 14-15kd

Preparation

This antibody was obtained by immunization of Galectin-7 recombinant protein (Accession Number: NM_001042507). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


Validations



mouse skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:110841(Galectin-7 antibody) at dilution of 1:1500



Immunohistochemical of paraffin-embedded human cervix using Catalog No:110841(Galectin-7 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human cervix using Catalog No:110841(Galectin-7

Galectin-7 antibody



Catalog Number: 110841

antibody) at dilution of 1:100 (under 40x lens)