

GALK1 antibody



Catalog Number: 110842

Product name

GALK1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GALK1 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, mouse heart tissue. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of GALK1 recombinant protein (Accession Number: NM_001381985). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110842(GALK1 antibody) at dilution of 1:1200