

# GALNT14 antibody



Catalog Number: 110846

## Product name

GALNT14 antibody

Recommended Dilution:

WB: 1:500-1:5000

## Specificity


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

## Antibody description

GALNT14 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, mouse kidney tissue, mouse lung tissue, rat testis tissue. Positive IHC detected in human kidney tissue, human heart tissue, human liver tissue, human ovary tissue, human skin tissue, human spleen tissue, human testis tissue. Observed molecular weight by Western-blot: 64kd

## Validations




## Preparation

This antibody was obtained by immunization of GALNT14 recombinant protein (Accession Number: NM\_024572). Purification method: Antigen affinity purified.

mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110846(GALNT14 antibody) at dilution of 1:400

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT


## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human kidney using Catalog No:110846(GALNT14 antibody) at dilution of 1:200 (under 10x lens)

## Ig Type

Rabbit IgG



## Applications

ELISA, WB, IHC

## Dilutions

# GALNT14 antibody



Catalog Number: 110846

---

Immunohistochemical of paraffin-embedded

human kidney using Catalog No:110846(GALNT14 antibody) at dilution of 1:200 (under 40x lens)