

GALNT4 antibody



Catalog Number: 110847

Product name

GALNT4 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

GALNT4 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IHC detected in human prostate cancer. Observed molecular weight by Western-blot: 67kd

Preparation

This antibody was obtained by immunization of GALNT4 recombinant protein (Accession Number: NM_003774). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110847(GALNT4 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:110847(GALNT4 antibody) at dilution of 1:100 (under 10x lens)