

# GALNTL5 antibody



Catalog Number: 110851

## Product name

GALNTL5 antibody

## Specificity

Human; other species not tested.

## Antibody description

GALNTL5 Rabbit Polyclonal antibody. Positive IHC detected in human lung tissue. Positive WB detected in human placenta tissue, human testis tissue. Observed molecular weight by Western-blot: 51 kDa

## Preparation

This antibody was obtained by immunization of GALNTL5 recombinant protein (Accession Number: XM\_006715862). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:110851(GALNTL5 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lung using Catalog No:110851(GALNTL5 antibody) at dilution of 1:100 (under 40x lens)