

GALP antibody



Catalog Number: 110852

Product name

GALP antibody

Specificity

Human; other species not tested.

Antibody description

GALP Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue, SH-SY5Y cells. Observed molecular weight by Western-blot: 5 kDa

Preparation

This antibody was obtained by immunization of GALP recombinant protein (Accession Number: NM_033106). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110852(GALP Antibody) at dilution of 1:300



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:110852(GALP Antibody) at dilution of 1:300