

gamma Actin antibody



Catalog Number: 110854

Product name

gamma Actin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

gamma Actin Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A549 cells, mouse kidney tissue, mouse lung tissue. Observed molecular weight by Western-blot: 42-45 kDa

Preparation

This antibody was obtained by immunization of gamma Actin recombinant protein (Accession Number: NM_001199954). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110854(ACTG1 antibody) at dilution of 1:500