

# gamma Actin antibody



Catalog Number: 110854

## Product name

gamma Actin antibody

Polyclonal

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat; other species not tested.

## Applications

ELISA, WB

## Antibody description

gamma Actin Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A549 cells, mouse kidney tissue, mouse lung tissue. Observed molecular weight by Western-blot: 42-45 kDa

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of gamma Actin recombinant protein (Accession Number: NM\_001199954). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110854(ACTG1 antibody) at dilution of 1:500

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality