

# GAN antibody

Catalog Number: 110861

---

**Product name**

GAN antibody

**Specificity**

human, mouse, rat; other species not tested.

**Antibody description**

GAN Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HEK-293 cells. Observed molecular weight by Western-blot: 68kd

**Preparation**

This antibody was obtained by immunization of GAN recombinant protein (Accession Number: NM\_022041). Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG


**Applications**

ELISA, WB

**Dilutions**

Recommended Dilution:

WB: 1:500-1:5000

**Validations**

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110861(GAN antibody) at dilution of 1:500