

GAS6 antibody



Catalog Number: 110872

Product name

GAS6 antibody

Specificity

Human; other species not tested.

Antibody description

GAS6 Rabbit Polyclonal antibody. Positive WB detected in DU 145 cells. Positive IHC detected in human lung tissue, human ovary tissue. Observed molecular weight by Western-blot: 75 kDa

Preparation

This antibody was obtained by immunization of GAS6 recombinant protein (Accession Number: NM_000820). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



DU 145 cells were subjected to SDS PAGE followed by western blot with Catalog No:110872(GAS6 Antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human lung using Catalog No:110872(GAS6 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human lung using Catalog No:110872(GAS6 antibody) at dilution of 1:50 (under 40x lens)