

# GATA3 antibody



Catalog Number: 110880

## Product name

GATA3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

GATA3 Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue, K-562 cells. Observed molecular weight by Western-blot: 45-50 kDa

## Preparation

This antibody was obtained by immunization of GATA3 recombinant protein (Accession Number: NM\_001002295). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:110880(GATA3 Antibody) at dilution of 1:600



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:110880(GATA3 Antibody) at dilution of 1:300