

# GATA4-Specific antibody



Catalog Number: 110882

## Product name

GATA4-Specific antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GATA4-Specific Rabbit Polyclonal antibody.  
Positive IHC detected in human testis tissue, human heart tissue. Positive IP detected in mouse heart tissue. Positive WB detected in mouse liver tissue, mouse heart tissue, rat heart tissue. Observed molecular weight by Western-blot: 45kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_002052).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000


## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110882(GATA4-Specific antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human testis using Catalog No:110882(GATA4-Specific antibody) at dilution of 1:50 (under 10x lens)



# GATA4-Specific antibody



Catalog Number: 110882

Immunohistochemical of paraffin-embedded human testis using Catalog No:110882(GATA4-Specific antibody) at dilution of 1:50 (under 40x

IP Result of anti-GATA4-Specific (IP:Catalog No:110882, 4ug; Detection:Catalog No:110882 1:600) with mouse heart tissue lysate 3200ug.

lens)

