

# GATA6 antibody



Catalog Number: 110884

## Product name

Polyclonal

GATA6 antibody

## Ig Type

## Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

GATA6 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, BxPC-3 cells, MCF7 cells. Observed molecular weight by Western-blot: 60-64kd

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_005257). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110884(GATA6 antibody) at dilution of 1:500