

# GBF1 antibody



Catalog Number: 110890

## Product name

GBF1 antibody

## Specificity

Human; other species not tested.

## Antibody description

GBF1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 206 kDa

## Preparation

This antibody was obtained by immunization of GBF1 recombinant protein (Accession Number: NM\_001391924). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110890(GBF1 Antibody) at dilution of 1:600



Immunofluorescent analysis of HeLa cells using Catalog No:110890(GBF1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)