

GCK antibody



Catalog Number: 110904

Product name

Polyclonal

GCK antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

GCK Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 52 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of GCK recombinant protein (Accession Number: NM_000162). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110904(GCK Antibody) at dilution of 1:600