

GCM2 antibody



Catalog Number: 110909

Product name

GCM2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GCM2 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, human brain tissue, mouse kidney tissue. Observed molecular weight by Western-blot: 50-56 kDa

Preparation

This antibody was obtained by immunization of GCM2 recombinant protein (Accession Number: NM_004752). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110909(GCM2 antibody) at dilution of 1:500