

# GCN5 antibody



Catalog Number: 110910

---

## Product name

GCN5 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

GCN5 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, mouse ovary tissue. Observed molecular weight by Western-blot: 94 kDa

## Preparation

This antibody was obtained by immunization of GCN5 recombinant protein (Accession Number: NM\_021078). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110910(KAT2A antibody) at dilution of 1:800