

GDAP1 antibody



Catalog Number: 110919

Product name

GDAP1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GDAP1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, fetal human brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 36-41 kDa

Preparation

This antibody was obtained by immunization of GDAP1 recombinant protein (Accession Number: NM_001040875). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110919(GDAP1 antibody) at dilution of 1:500