

GDE1 antibody



Catalog Number: 110922

Product name

GDE1 antibody

Immunogen

[Human GDE1 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GDE1 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, HEK-293 cells. Observed molecular weight by Western-blot: 55-66kd

Preparation

This antibody was obtained by immunization of GDE1 recombinant protein (Accession Number: NM_016641). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110922(GDE1 antibody) at dilution of 1:500