

GDI1 antibody



Catalog Number: 110926

Product name

Polyclonal

GDI1 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

GDI1 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, human brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 51-55kd

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001493). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:110926(GDI1 antibody) at dilution of 1:1000