

GDI1 antibody



Catalog Number: 110926

Product name

GDI1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

GDI1 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, human brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 51-55kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001493). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:110926(GDI1 antibody) at dilution of 1:1000