

GDPD2 antibody



Catalog Number: 110929

Product name

GDPD2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GDPD2 Rabbit Polyclonal antibody. Positive WB detected in mouse skin tissue. Observed molecular weight by Western-blot: 62 kDa

Preparation

This antibody was obtained by immunization of GDPD2 recombinant protein (Accession Number: NM_017711). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:110929(GDPD2 Antibody) at dilution of 1:600