

GFRA2 antibody

Catalog Number: 110954

Product name

GFRA2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GFRA2 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, mouse brain tissue. Positive IHC detected in human liver tissue. Observed molecular weight by Western-blot: 34-40kd

Preparation

This antibody was obtained by immunization of GFRA2 recombinant protein (Accession Number: XM_011544484). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



WB result of Catalog No:110954(GFRA2




HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110954(GFRA2 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human liver using Catalog No:110954(GFRA2 antibody) at dilution of 1:50 (under 10x lens)

GFRA2 antibody

Catalog Number: 110954



Immunohistochemical of paraffin-embedded human liver using Catalog No:110954(GFRA2 antibody) at dilution of 1:50 (under 40x lens)