

GFRA3 antibody



Catalog Number: 110955

Product name

Polyclonal

GFRA3 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

GFRA3 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: 45 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of GFRA3 recombinant protein (Accession Number: NM_001496). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:110955(GFRA3 Antibody) at dilution of 1:600