

GGA2 antibody



Catalog Number: 110957

Product name

GGA2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GGA2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human kidney tissue, human liver tissue. Observed molecular weight by Western-blot: 67 kDa

Preparation

This antibody was obtained by immunization of GGA2 recombinant protein (Accession Number: NM_015044). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110957(GGA2 antibody) at dilution of 1:1000