

GGT5 antibody

Catalog Number: 110961

Product name

GGT5 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

GGT5 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 80kd

Preparation

This antibody was obtained by immunization of GGT5 recombinant protein (Accession Number: NM_004121). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110961(GGT5 antibody) at dilution of 1:400