

GH2 antibody

Catalog Number: 110963

Product name

GH2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GH2 Rabbit Polyclonal antibody. Positive IHC detected in human placenta tissue, human pituitary adenoma tissue. Positive WB detected in human placenta tissue. Observed molecular weight by Western-blot: 23-25kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002059). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:110963(GH2 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human placenta using Catalog No:110963(GH2 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human placenta using Catalog No:110963(GH2 antibody) at dilution of 1:50 (under 40x lens)