

Product name

GHR antibody

Specificity

Human, Pig; other species not tested.

Antibody description

GHR Rabbit Polyclonal antibody. Positive WB detected in pig lung tissue. Observed molecular weight by Western-blot: 120 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000163). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations

pig lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:110966(GHR Antibody) at dilution of 1:300