

# GINS3 antibody



Catalog Number: 110976

## Product name

GINS3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

GINS3 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 24kd

## Preparation

This antibody was obtained by immunization of GINS3 recombinant protein (Accession Number: NM\_022770). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110976(GINS3 antibody) at dilution of 1:800