

# GJB6 antibody

Catalog Number: 110985

## Product name

GJB6 antibody

## Immunogen

[Human GJB6 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

GJB6 Rabbit Polyclonal antibody. Positive IHC detected in human oesophagus tissue, human skin cancer tissue.

## Preparation

This antibody was obtained by immunization of GJB6 recombinant protein (Accession Number: NM\_001110221). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, IHC


## Dilutions

Recommended Dilution:  
IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human oesophagus tissue slide using Catalog No:110985(GJB6 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human oesophagus tissue slide using Catalog No:110985(GJB6 Antibody) at dilution of 1:50 (under 40x lens)