

GNAQ antibody

Catalog Number: 110991

Product name

GNAQ antibody

Specificity

human, mouse; other species not tested.

Antibody description

GNAQ Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human brain tissue, Jurkat cells, mouse brain tissue, mouse lung tissue, mouse skin tissue. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of GNAQ recombinant protein (Accession Number: NM_002072). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110991(GNAQ antibody) at dilution of 1:1000