

GNAQ antibody



Catalog Number: 110991

Product name

Polyclonal

GNAQ antibody

Ig Type

Specificity

Rabbit IgG

human, mouse; other species not tested.

Applications

Antibody description

ELISA, WB

GNAQ Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human brain tissue, Jurkat cells, mouse brain tissue, mouse lung tissue, mouse skin tissue. Observed molecular weight by Western-blot: 42 kDa

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Preparation

Validations

This antibody was obtained by immunization of GNAQ recombinant protein (Accession Number: NM_002072). Purification method: Antigen affinity purified.


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110991(GNAQ antibody) at dilution of 1:1000