

# GNB3 antibody



Catalog Number: 110997

## Product name

GNB3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GNB3 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 35-37 kDa

## Preparation

This antibody was obtained by immunization of GNB3 recombinant protein (Accession Number: NM\_002075). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110997(GNB3 antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells using Catalog No:110997(GNB3 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)