

GNG2 antibody



Catalog Number: 111001

Product name

GNG2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GNG2 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, mouse brain tissue. Observed molecular weight by Western-blot: 8 kDa

Preparation

This antibody was obtained by immunization of GNG2 recombinant protein (Accession Number: XM_017021377). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:111001(GNG2 antibody) at dilution of 1:400