

## Product name

GNG4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GNG4 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 8 kDa

## Preparation

This antibody was obtained by immunization of GNG4 recombinant protein (Accession Number: NM\_001098722). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111002(GNG4 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:111002(GNG4 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:111002(GNG4 antibody) at dilution of 1:100

# GNG4 antibody



Catalog Number: 111002

---

(under 40x lens)