

GNPNAT1 antibody



Catalog Number: 111010

Product name

GNPNAT1 antibody

Immunogen

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GNPNAT1 Rabbit Polyclonal antibody. Positive IP detected in mouse large intestine tissue. Positive WB detected in mouse large intestine tissue. Positive IHC detected in human liver tissue. Observed molecular weight by Western-blot: 21-23 kDa

Preparation

This antibody was obtained by immunization of GNPNAT1 recombinant protein (Accession Number: NM_198066). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions

Recommended Dilution:

WB: 1:200-1:1000


IP: 1:200-1:1000

IHC: 1:20-1:200

Validations




mouse large intestine tissue were subjected to SDS PAGE followed by western blot with Catalog No:111010(GNPNAT1 antibody) at dilution of 1:300




Immunohistochemical of paraffin-embedded human liver using Catalog No:111010(GNPNAT1 antibody) at dilution of 1:100 (under 10x lens)

GNPNAT1 antibody

Catalog Number: 111010



Immunohistochemical of paraffin-embedded human liver using Catalog No:111010(GNPAT1 antibody) at dilution of 1:100 (under 40x lens)



IP Result of anti-GNPAT1 (IP:Catalog No:111010, 3ug; Detection:Catalog No:111010 1:300) with mouse large intestine tissue lysate 3200ug.