

# GLDC antibody

Catalog Number: 111020

## Product name

GLDC antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

GLDC Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 113 kDa

## Preparation

This antibody was obtained by immunization of GLDC recombinant protein (Accession Number: NM\_000170). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111020(GLDC Antibody) at dilution of 1:1000



IP Result of anti-GLDC (IP:Catalog No:111020, 4ug; Detection:Catalog No:111020 1:1000) with mouse liver tissue lysate 4000ug.