

GLRA2 antibody



Catalog Number: 111031

Product name

GLRA2 antibody

WB: 1:300-1:1200

IHC: 1:20-1:200

Immunogen

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GLRA2 Rabbit Polyclonal antibody. Positive WB detected in PBMC cells. Positive IHC detected in mouse brain tissue. Observed molecular weight by Western-blot: 48 kDa

Preparation

This antibody was obtained by immunization of GLRA2 recombinant protein (Accession Number: NM_001118885). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

Validations




WB result of GLRA2 antibody (Catalog No:111031, 1:600) with PBMC cells. Proteintech α GLRA2 Catalog No:111031 in 10% BSA / TBS 0.05% Tween-20 O/N @ 4°C Goat anti-rabbit HRP 1:10000 in 5% milk / TBS 0.05% Tween-20 1h @



RT.

Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:111031(GLRA2 Antibody) at dilution of 1:50 (under 10x lens)



GLRA2 antibody



Catalog Number: 111031

Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog

No:111031(GLRA2 Antibody) at dilution of 1:50 (under 40x lens)