

# GLRA2 antibody



Catalog Number: 111031

## Product name

GLRA2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GLRA2 Rabbit Polyclonal antibody. Positive WB detected in PBMC cells. Positive IHC detected in mouse brain tissue. Observed molecular weight by Western-blot: 48 kDa

## Preparation

This antibody was obtained by immunization of GLRA2 recombinant protein (Accession Number: NM\_001118885). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:300-1:1200

IHC: 1:20-1:200


## Validations



WB result of GLRA2 antibody (Catalog No:111031, 1:600) with PBMC cells. Proteintech  $\alpha$ GLRA2 Catalog No:111031 in 10% BSA / TBS 0.05% Tween-20 O/N @ 4°C Goat anti-rabbit HRP 1:10000 in 5% milk / TBS 0.05% Tween-20 1h @



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:111031 (GLRA2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog

# GLRA2 antibody



Catalog Number: 111031

---

No:111031(GLRA2 Antibody) at dilution of 1:50 (under 40x lens)