

GLRB antibody



Catalog Number: 111033

Product name

GLRB antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GLRB Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, mouse spinal cord tissue. Observed molecular weight by Western-blot: 56 kDa

Preparation

This antibody was obtained by immunization of GLRB recombinant protein (Accession Number: NM_001166060). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111033(GLRB Antibody) at dilution of 1:1000