

# GOLGA3 antibody



Catalog Number: 111043

## Product name

GOLGA3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

GOLGA3 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, HEK-293 cells, HeLa cells, human heart tissue, human testis tissue, PC-3 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 167kd

## Preparation

This antibody was obtained by immunization of GOLGA3 recombinant protein (Accession Number: NM\_001172557). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

## Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:111043(GOLGA3 antibody) at dilution of 1:2000



Immunofluorescent analysis of HepG2 cells, using GOLGA3 antibody Catalog No:111043 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).