

GOSR2/Membrin antibody



Catalog Number: 111049

Product name

GOSR2/Membrin antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GOSR2/Membrin Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in Jurkat cells, human liver tissue, mouse liver tissue, mouse stomach tissue. Observed molecular weight by Western-blot: 25 kDa

Preparation

This antibody was obtained by immunization of GOSR2/Membrin recombinant protein (Accession Number: NM_004287). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111049(GOSR2 antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells using Catalog No:111049(GOSR2 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG