

GPI antibody

Catalog Number: 111064

Product name

GPI antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GPI Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, HeLa cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 63 kDa

Preparation

This antibody was obtained by immunization of GPI recombinant protein (Accession Number: NM_000175). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:111064(GPI antibody) at dilution of 1:500



IP Result of anti-GPI (IP:Catalog No:111064, 4ug; Detection:Catalog No:111064 1:500) with HeLa cells lysate 1200ug.