

# GPR103 antibody

Catalog Number: 111070

## Product name

GPR103 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

GPR103 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 50 kDa

## Preparation

This antibody was obtained by immunization of GPR103 recombinant protein (Accession Number: BC128133). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:50-1:500


## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111070(QRFPR Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:111070(QRFPR Antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:111070(QRFPR Antibody) at dilution of 1:100 (under 40x lens).