

GLYAT antibody



Catalog Number: 111082

Product name

Polyclonal

GLYAT antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

GLYAT Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 35 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of GLYAT recombinant protein (Accession Number: NM_005838). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111082(GLYAT Antibody) at dilution of 1:600